

Note: The most current version of these instructions can be found online at www.grahamfield.com

IMPORTANT SAFETY INFORMATION: PLEASE READ

- ⚠ **WARNING: Important! Read and understand these instructions before installing or using the GF13702 Advanced Patient Alarm. If you do not understand any part of these instructions, contact a healthcare professional for direction in the use of this product. If the patient alarm is not properly set up and adjusted, personal injury and/or damage to the patient alarm could result.**
- ⚠ **WARNING: If components are damaged or missing, contact your dealer immediately. DO NOT use substitute parts. Use only Lumex® replacement parts. Non-Lumex® replacement parts could cause personal injury and/or damage to patient alarm.**
- ⚠ **WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper assembly or use of this product.**

UNPACKING

Check for any obvious damage to the carton or its contents. If damage is evident, DO NOT USE. Contact carrier / distributor for further instruction.

INTENDED USE


The intended use of the GF13702 Advanced Patient Alarm is to alert caregivers when a patient attempts to exit a wheelchair, chair or bed without assistance. Removal of pressure from the sensor pad beneath the patient activates the alarm.

Contraindications

- ⚠ **WARNING: This device is not a substitute for routine visual monitoring by caregiver, and is intended to be used in conjunction with a comprehensive fall management program.**
- ⚠ **WARNING: No claim is made that this device will prevent patient falling or reduce the likelihood of injury as a result. The caregiver is responsible for ensuring that this product is properly installed, tested, and functioning before each use.**
- ▲ **CAUTION: This device is not suitable for all patients, including but not limited to, those wishing to defeat the alarm system.**
- ▲ **CAUTION: This product is designed for individuals weighing 74 lb or more.**

- ▲ **CAUTION: Seat cushions, sling seats, and other environmental factors may negatively affect alarm operation.**

ALARM SETUP



1. Remove the screw that holds the battery compartment cover in place.
2. Inside the battery compartment (see Fig. B):
 - a. Set the volume switch to LO, HI, or RE (remote) volume. For description of the remote alarm feature, see next page.
 - b. Set the delay switch to either 0 or 2 seconds. Selecting 2-second alarm delay may reduce the likelihood of false alarms when used with restless patients.
 - c. Install a new 9V alkaline battery (not included).
 - d. Replace battery compartment cover and reinstall screw.

Alarm Testing

- ⚠ **WARNING: Test alarm before each use.**
1. Insert the sensor pad's telephone-style plug into the monitor's sensor mat jack (see Fig. C) until it clicks into place.
 2. To protect the cable plug end from accidental pull force, route the cable through the strain relief (see Fig. C).
 3. Use your hand to apply pressure to the sensor pad. The monitor will automatically sense the pressure and beep when activated.
 4. Remove the pressure and the alarm will sound. To silence the alarm, press the reset button on the side of the monitor (see Fig. D).

Patient Setup

1. Follow this procedure before each use.
 - a. If using a chair sensor pad, center the chair sensor pad under the patient's buttocks.
 - b. If using a bed sensor pad, center the bed sensor pad under the patient's buttocks.

2. The monitor will automatically sense the pressure and beep when activated.
3. Have the patient attempt to stand to remove pressure from the sensor pad. The alarm will sound. To silence the alarm, press the reset button on the side of the monitor (see Fig. D).
4. **Place monitor out of reach of patient.**

⚠ WARNING: If alarm fails to sound, immediately discontinue use and review alarm setup.

Note: False alarms are often the result of poor pad placement. For instance, patients may position themselves differently on a bed with a raised head section than on a completely flat mattress.

To reduce the likelihood of false alarms, the caregiver *must* ensure that the pad is correctly placed for each patient's unique positioning needs.

Remote Alarm (Optional)

The remote alarm option allows the caregiver to silence the local alarm and sound the alarm remotely through the nurse call system.

1. Select the remote alarm option: Open the battery compartment and move the volume switch to "RE" (remote) position.
2. Connect the monitor to the nurse call system using an adapter cable and splitter (not included).

⚠ WARNING: Test the alarm in the nurse call system before each use. If alarm fails to sound, immediately discontinue use and review alarm setup.

Note: Optional accessories (adapter cable and AC adapter) are available through GF Health Products, Inc.

MAINTENANCE

⚠ WARNING: Test the alarm before each use. If alarm fails to sound, remove it from service immediately.

Check before each use that all alarm components are functional and battery is charged. To replace battery, see ALARM SETUP section on previous page.

Note: The low battery indicator on the front of the alarm will flash when battery is low and needs to be replaced. Replace battery immediately if indicator flashes (see Fig. A).

Cleaning

1. To clean: Gently wipe alarm components with a soft, clean cloth.

▲ CAUTION: Do not use cleansers that could damage the alarm's finish. Do not submerge any alarm components in liquid.

WARRANTY

GF Health Products, Inc. ("Graham-Field") offers a one-year limited warranty against manufacturer's defects on the Lumex® GF13702 Advanced Patient Alarm.

The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alterations or disassembly by unauthorized individuals.

If a product is deemed to be under warranty, GF Health Products, Inc. shall provide, at its option, (1) replacement of any defective part or product or (2) a credit of the original selling price made to GF Health Products, Inc.'s initial customer. The warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to GF Health Products, Inc.

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

Manufactured for
GF Health Products, Inc.
Atlanta, GA 30360



GRAHAM-FIELD
www.grahamfield.com

GF Health Products, Inc. is not responsible for typographical errors. Packaging, warranties, products and specifications are subject to change without notice. Graham-Field and Lumex are registered trademarks of GF Health Products, Inc.

Manufactured by Arrowhead Healthcare Supply, Duluth, MN 55802 • Made in China

© 2010 GF Health Products, Inc. • 2935 Northeast Parkway • Atlanta, GA 30360
tel 800-347-5678 • fax 800-721-0601